

The `visrock` user manual

Daniel Lübbert

July 23, 2006

(Id: visrockManual.tex,v 1.1 2005/08/02 14:54:41 luebbert Exp)

Contents

1	Introduction - what the <code>visrock</code> program can be used for	3
1.1	Structure of this manual	5
2	Preconditions to using <code>visrock</code>	6
2.1	Obtaining <code>visrock</code>	6
2.2	License policy	6
2.3	System requirements	6
2.4	Installation	7
2.5	Launching <code>visrock</code>	7
3	Using <code>visrock</code>, I: Basic functionality – for beginners	9
3.1	Data input (reading images)	9
3.2	Image analysis options	11
3.3	Saving the results to disk	14
3.4	Terminating the program	15
4	Using <code>visrock</code>, II: Advanced functionality – for experienced users	16
4.1	Curve analysis options	16
4.2	Curve plotting options	16
4.3	Image Plotting options	18
4.4	Axis units	18
4.5	Histogram display	22
4.6	Line profile extraction from images	23
4.7	Correlations	23
4.8	Noise reduction	24
4.9	Image processing options - sharpening, smoothing and peak detection	29

5	Using visrock, III: Special functions for specific applications	32
5.1	Functions for Rocking Curve Imaging	32
5.2	Functions for Tomography	33
5.3	Functions for Fine-Slicing	34
5.4	Functions for Small-angle scattering and general Crystallography	35
5.5	Functions for Detector Calibration	35
6	Planned future developments of visrock	38
7	List of known bugs	39
7.1	Bugs	39
7.2	Not bugs - features	39
8	Impressum	41
9	Literature references	42

1 Introduction - what the visrock program can be used for

The `visrock` program was written to facilitate the interactive handling and analysis of digital X-ray images. It was originally developed for applications in Rocking Curve Imaging, but was later extended to several other fields of X-ray imaging. It can be useful in various applications that are based on linear sequences of 2D images ("sequential imaging"):

Rocking Curve Imaging (RCI): Visualizing the spatial distribution of mosaicity and lattice quality by a series of plane-wave X-ray diffraction topographs recorded along the rocking curve of a single crystal [LBH00, MLK03]

X-ray tomography: 3D image reconstruction from a series of projection radiographs [Her80]. `visrock` cannot do the reconstruction, but can perform a few tasks for preparing experimental data for reconstruction and for checking data integrity.

Fine-Slicing (or "Super-fine Φ slicing"): Measuring the mosaicity of (mostly protein) crystals along all directions of reciprocal space by simultaneously recording rocking curves from a great number of crystallographic reflections, using a CCD camera covering a wide range of exit angles [BSL00].

Diffraction-enhanced imaging (DEI): Recording projection radiographs of a (possibly non-crystalline) sample with a crystal analyzer as an additional optical element in the exit beam, at a series of angular positions along the rocking curve of the analyzer crystal (see e.g. [Pag03]).

Reticulography: Visualizing lattice irregularities in a crystal by monitoring the distortion of a rectangular grid pattern in X-ray diffraction topographs recorded at several different sample-camera distances [LM96].

Detector calibration: `visrock` includes features for testing and quantifying the characteristics of CCD cameras (and other digital area detectors) in terms of reproducibility and noise level, linearity, sensitivity, and spatial homogeneity of response. These can be determined from various series of dark field and flat field images.

The `visrock` program is a graphical user interface (GUI) – it is entirely based on a point-and-click strategy. The main aim is to simplify tasks that occur frequently and repetitively in image sequence analysis. No keyboard input of any kind is required – a simple mouse will suffice as an input device. Among the functions that can be executed in this way are (sorted roughly from simple to more complex):

- Reading image data from a large variety of file formats.
- Color plotting of (scalar, i.e. intrinsically "black-and-white") images to screen or to postscript output files

- Interactively adjustable color scales, either linear or logarithmic
- Plot range selection based on absolute numbers or on percentile values
- Playing movies of image sequences
- Format conversion of grayscale images (e.g. from *.edf to *.png, from *.tif to *.edf, and many more) – for single images or entire image series
- Interactive zooming (specifying areas of interest with the mouse)
- Monitoring the evolution of the intensity value in individual pixels along an image sequence (file series), or of the average intensity of a group of pixels (Interactive extraction of local or regional "Rocking Curves")
- Automatic peak detection in extracted profiles
- Determination of peak half widths (FWHM) in profiles
- Gaussian fitting of profile peaks
- Line profile extraction from images, in horizontal or vertical direction
- Horizontal and vertical projections (summation) of image contents
- Image rotations, by arbitrary angles
- Image statistics
- Dark field and flat field corrections
- Renormalization of an image sequence, based on the intensity statistics in a small region of interest/ROI (e.g. to compensate for a decay in primary intensity)
- Rebinning (in 1, 2 or 3 dimensions)
- Sharpening and smoothing of images
- Noise reduction
- 2D peak detection in images
- Fourier transforms of images (in 2 dimensions)

1.1 Structure of this manual

The remainder of this manual is divided into 4 main chapters:

Chapter 2 describes how to prepare for the first use of **visrock**.

Chapter 3 discusses a first complete round of using **visrock**, from data input to result output. It is limited to the absolute essentials, and meant to be short and easy to read. It should be considered as obligatory reading by every new user of the program.

Chapter 4 provides guidance to the more advanced features of the **visrock** program. It covers in detail all the available image analysis, processing and plotting options, and serves as a reference to those users who want to get the most out of the program.

Finally, chapter 5 covers a few very specific functions related to single scientific applications. It shows how individual task occurring repetitively in some fields of X-ray imaging may be carried out by **visrock**.

2 Preconditions to using visrock

2.1 Obtaining visrock

The `visrock` software can be obtained directly from the author, based on a personal collaboration agreement.

2.2 License policy

The `visrock` software may be used free of cost and for an unlimited time. However, it must not be redistributed, nor made publicly available on the internet, on an intranet or on any other publicly available storage space.

The author will be happy to receive suggestions for new collaborators and users of the `visrock` project.

2.3 System requirements

2.3.1 Software requirements

`visrock` was written in the IDL programming language (see <http://www.rsinc.com>). The software is distributed as an IDL `*.sav` file, which requires a working IDL installation in order to be run. Ideally, the user's institution should have an IDL license – which is the case at almost all synchrotrons, many universities and other research institutions. However, from IDL version 6.0 on, `*.sav` files can also be run from unlicensed IDL installations (the so-called "IDL Virtual Machine").

By default, `visrock` is optimized for the latest version of IDL (version 6.0, as of July 2004).

2.3.2 Operating system

IDL runs on almost any modern operating system (Linux, Unix'es, all kinds of Windows, Apple, ...). `visrock` can therefore be run on any type of computer. Only some special conveniences realized by calls to external programs are limited to Linux/UNIX systems. The concerns in particular the input of compressed data and the display on screen of postscript files.

2.3.3 Hardware requirements

RAM: `visrock` assumes that the entire sequence of images can be loaded into RAM simultaneously. This is sometimes a very severe restriction to the allowed number of images and/or pixels.

An example: Even for a well-equipped computer with 2 GB RAM and for images of 2048^2 pixels à 2 Bytes (typical case for ESRF-Frelon cameras), this means that

no more than 200–230 images can be analyzed at once (leaving some RAM for IDL itself and other system applications).

Larger sequences can either be analyzed piecewise (e.g. in four "quadrants"), or on systems with even more memory (such as the ESRF-*indigo* machines with 16 GB of RAM).

Screen: `visrock` can be run on screens having at least 1024×768 pixels, but is optimized for larger screens (1600×1200 is ideal).

2.4 Installation

`visrock` does not need to be installed in any way. Just copy the downloaded file `visrock.sav` to the directory of your choice - that's it.

2.5 Launching `visrock`

2.5.1 With an IDL license

Start IDL by typing "idl" on the Unix prompt, or by clicking on the appropriate icon under Windows.

On the IDL command line, type

```
restore, '<path>/visrock.sav'
```

where *< path >* is the path to the directory where you saved `visrock.sav`.

Provided that you have correctly set the `IDL_PATH` environment variable, you should now see the `visrock` GUI (see Fig. 1).

2.5.2 Without an IDL license

Starting from IDL version 6.0 on, precompiled IDL programs (such as `visrock.sav`) may be run in the "IDL Virtual Machine (VM)" even without an IDL license. For further explanations, see the file `idlv.m.pdf` provided in the `help` subdirectory of your IDL distribution.

To start `visrock` in the IDL VM, proceed as follows:

Under Linux: On your command line, type

```
idl -vm=<path>/visrock.sav
```

where *< path >* is the directory where you stored `visrock.sav`.




Figure 1: Appearance of the visrock GUI, in Linux style. The two plot areas are for image display (left) and curve display (right). The meaning of the remaining buttons and fields will be explained in subsections 3 and 4.

Under Windows: Drag-and-drop the `visrock.sav` program onto the IDL-VM icon on your desktop.

In both cases, you will have to click away the IDL-VM splash screen.

Please be aware that a few limitations exist in the unlicensed mode of IDL, as described in `idlv.m.pdf`. In particular, IDL's `execute` function will not work, which means that commands entered via the "execute" text field at the bottom of the `visrock` window cannot be carried out.

3 Using visrock, I: Basic functionality – for beginners

Basic usage of `visrock` usually includes three steps: Data input, image analysis, and output of results (most often as postscript plots).

3.1 Data input (reading images)

The `visrock` program works on a 3D data set (intensity data as function of e.g. x , y and ω). These are usually read into memory as a sequence of 2D images (intensity data as a function of x and y , for one value of ω at a time).

3.1.1 Reading individual images

Single images can be read from hard disk by clicking on the `Read CCD-files` button, in the upper right of the `visrock` GUI. Multiple images can be selected by using the `SHIFT` and/or `CTRL` keys, and will be read in sequence after clicking `OK`.

3.1.2 Reading large series of images

When reading large series of images, it may be more convenient to specify the images to be read via a `*.txt` file containing the names of all the image files, in a single column. This can be done by clicking on the `Read File list` button.

3.1.3 Supported image formats

The format of the image files will be automatically analyzed by `visrock` from the file name extension. The following image formats are understood by `visrock`:

- ***.tiff**: Standard digital image format.
- ***.jpg**: Standard digital image format.
- ***.gif**: Standard digital image format.
- ***.png**: Standard digital image format.
- ***.edf**: ESRF Data Format, produced by 2D detectors such as the "Frelon" cameras. ASCII header, plus uncompressed binary image data.
- ***.mccd**: Image format produced by MarResearch CCD cameras, in widespread use in protein crystallography. Regular tiff format, plus extra binary header.
- ***.b16**: Slightly exotic digital image format, apparently produced by CCD cameras from PCO.
- ***.gqe**: Binary image format produced by the Spectra program, developed by Thorsten Kracht and used at many HASYLAB beamlines.

***.xdr:** Even more exotic general file format, somehow device-independent, read and written by older IDL programs (in use at BW4@HASYLAB).

In addition to the above (and as a default, if the file name matches the standard for none of the above), **visrock** can analyze a general binary image format, consisting of a H-Byte header (which will be ignored) plus a rectangular array of $N \times M$ pixels, with L Bytes per pixel. The integer numbers H, N, M and L will be requested interactively until the expected number of Bytes ($N * M * L + H$) matches the actual file size.

As an extra feature, **visrock** can read compressed (zipped) copies of all available image formats. This will work primarily on Linux systems; it may or may not work under Windows. It is particularly useful with ***.edf** images, which are uncompressed in their native form (unlike ***.jpg** or ***.png**) and therefore waste enormous amounts of hard disk space.

Currently, ***.gz** and ***.bz2** compressed formats are supported. When **visrock** encounters a compressed image, it is first uncompressed to a temporary file, and then read as usual. Decompression is usually rather fast (some tenths of a second), even if the original compression took considerably longer.

Note that for this to work, two conditions must be fulfilled:

- The appropriate decompression program must be found in the PATH (**gunzip** or **bunzip2**).
- You must have write permission on the **/tmp** directory.

3.1.4 Input details

After selection of the files to be read, a few more details of the reading process can be specified in a separate window (see Fig. 2). From top to bottom, the most important choices are (all parameters must be integer numbers):

bin_factor: Rebin original images, i.e. reduce their number of pixels via averaging, by this factor. The *x*- and *y*- dimension will be treated alike, i.e. reduced by the same factor. Must be an integer.

use every Nth file only: Read only image files nr. 0, *N*, 2*N*,... and disregard files nr. 1 – (*N* – 1) etc. Helps to save RAM, and considerably speeds up the reading process (by *at least* a factor *N*). Must be an integer.

wait time: If the next file to be read is not found on hard disk (e.g. because it is just being read out from the camera), wait *x* seconds and try again before giving up. Helpful when starting a reading process while the series of images is still being recorded.

x-start, x-stop, y-start, y-stop: Keep only a restricted area of the original images in memory, and discard the outer regions. Helpful in order to save memory space when the features of interest occupy only parts of the image surface. Values must be $\in [0, N - 1]$.




Figure 2: Popup window allowing to specify details of the image reading process.

The remaining fields of the window shown in Fig. 2 will be discussed in Chapter 5.

3.2 Image analysis options

In the left drawing area, 2D maps can be plotted which represent sections or projections from the 3D data block that resides in the program's memory. A large number of options is available, as can be seen from the list below. The first three options represent sections from the original data, whereas all the remaining options represent composite images generated from the entire stack of data, mostly by projection along the third (angular) dimension.

x-y (Single images): Show one of the original images. Its number is selectable via the slider at the bottom of the GUI.

x-z (x-omega): Show a section from the 3D data block, perpendicular to the original image planes. This is a sinogram in the case of tomography with vertical rotation axis (along y).

y-z (y-omega): Show a section from the 3D data block, perpendicular to the original image planes. This is a sinogram in the case of tomography with vertical rotation axis (along x).

integrated Intensity: For each pixel, show the sum of intensities across the entire image stack.

maximum Intensity: For each pixel, show the maximum intensity reached anywhere across the entire image series.

maxpos: For each pixel, show the angular position (or image number) where the maximum is reached.

minimum Intensity: For each pixel, show the minimum of intensity reached anywhere across the image series. Rarely used; sometimes useful with tomography images.

minpos: For each pixel, show the angular position (or image number) where the minimum is reached.

Int_and_Pos: Under development, no documentation yet. Don't use.

FWHM: For each pixel, calculate the full-width-at-half-maximum of the main peak in the local intensity curve.

FW_percent_M: For each pixel, calculate the full width at $x\%$ of peak height, with selectable percentage x . For $x = 50$, the result is identical to **FWHM**. Sometimes more representative than the **FWHM** map in case of multiple peaks or peaks with broad tails.

ratio FWHM / FW_percent_M: Quotient of the latter two quantities. Sometimes useful to detect regions with different numbers of peaks (overlapping grains in RCI).

Lintegr.d-restricted range: Summation of intensity across a restricted range of image numbers only, rather than the entire series. Will ask for numbers of first and last image. Useful when multiple peaks exist for each pixel.

maximum-restricted range: Similarly for maximum intensity.

maxpos-restricted range: Similarly for position of maximum.

FWHM-restricted range: Similarly for peak FWHM.

Angular center of mass: For each pixel, calculate the center of mass of the associated intensity distribution as a function of angle or image number (the $I(\omega)$ distribution). Result is in units of angle. Alternative to **maxpos**, statistically more robust: The resulting angle is slightly less sensitive to noise than **maxpos**.

Angular standard deviation: For each pixel, calculate the standard deviation of the local $I(\omega)$ distribution. Result is in units of angles, and is a measure of peak width. It is a statistically more robust alternative to **FWHM**, which may be preferable in case of noisy data, multiple peaks, and/or peaks with broad tails.

N_peaks: Calculate the number of peaks in each local intensity curve. Calculation based on the algorithm described in [Sil96].

Average: For each pixel, calculate the average intensity across the image series.

Median: For each pixel, calculate the median intensity across the image series.

max-min: For each pixel, calculate the difference between maximum and minimum intensity across the image series. Sometimes useful for tomography data.

Standard deviation: For each pixel, calculate the standard deviation of the associated vector of intensity values. Result is in terms of intensities. High values indicate that something changes somewhere in the entire series.

I over sigma: Ratio of average intensity over intensity standard deviation for each pixel.

noise <=> Poissonian?: Intensity standard deviation divided by square root of average intensity, for each pixel. Useful when evaluating the homogeneity of a series of flat field images.

wafer_present_yes_no: Binary image (yes/no) indicating which pixels “looked” on diffracting material (yes) or air (no). Calculation based on selectable upper/lower limits for **FWHM** and **max_Int**. Must have been calculated before using the **mask wafer area** button, located on the right side of the **visrock** GUI.

threshold: Binary image (yes/no) showing where **maximum Intensity** exceeded a selectable threshold.

Gauss_fit max_Int: Apply Gaussian fit to (highest) peak in each pixel’s intensity distribution, and display the peak height of the fitted curve. Robust alternative to **max_Int**, less sensitive to noise.

Gauss_fit max_Pos: Similarly, display angular peak position of fitted curve.

Gauss_fit FWHM: Similarly, display half width of fitted curve. Very interesting alternative to **FWHM**: less sensitive to noise, *and* less prone to errors due to discrete scan step width.

Linear fit: Background: Apply linear fit to local intensity distribution, and display residual (fitted intensity for $\omega = 0$). Useful mainly for detector calibration purposes.

Linear fit: Slope: Similarly, display slope of linear curve.

Linear fit: Chi²: Display chi-square for local fit. High values indicate deviations from linear local curve shape.

Linear fit: Status: Check result of local fit (0 for success, 1 for error, 2 for inaccurate result).

3D-Intensity: Display the entire 3D data block, similarly to volume data in tomography.

3D-Peak locations: Undocumented feature.

Principal components: Apply principal components analysis to the image stack. This is a standard procedure in multivariate (multi-dimensional) statistics: A new, rotated coordinate system is introduced in such a way that the first components of the (new) data vectors contain a maximum of the total variance of the data, and the last component a minimum.

For application to an image series, the image number is regarded as the “component index” of the data “vector”, and each pixel as an independent measure of the same “vectorial” quantity. The “coordinate transformation” creates a new set of images, the first of which has a maximum of contrast, whereas the last ones are rather “flat”.

This algorithm is a rather efficient method for contrast enhancement, but somehow dissatisfactory for *quantitative* imaging in the sense that the intensity values of the resulting images do not have any clearly defined unit.

3.3 Saving the results to disk

3.3.1 Postscript output

As soon as satisfactory plots of the relevant quantities have been obtained in the two plot areas in the **visrock** GUI, the images can be output almost identically to postscript files. To that end, click the **Map- > *.ps** and/or **Curve- > *.ps** buttons in the top row of the **visrock** GUI. A default filename will be suggested, containing quite detailed information on the kind of plot. The type of postscript plot (encapsulated yes/no) will be adapted to the final filename chosen (*.ps or *.eps).

On appropriate operating systems (Linux/UNIX), the newly created postscript file will be opened in **ghostview** for immediate checking/printing.

3.3.2 Data output

The results displayed in the plotting areas may as well be saved as raw data, for later display with other programs than **visrock**. To that end, click the **Map- > *.img** and/or **Curve- > *.dat** buttons in the top row of the **visrock** GUI. Maps can be stored in (almost) any of the image formats discussed in section 3.1.3. Curves will be saved as ascii files, with two data columns and one row for each image in the series.

3.4 Terminating the program

When your work is done, terminate the program by clicking either the **Quit** or the **Exit** button in the upper left of the **visrock** GUI. The former will close the GUI, but leave the IDL session open, with the values of various variables still in memory. The latter will close the IDL program as well and exit to the Linux shell / Windows desktop.

4 Using visrock, II: Advanced functionality – for experienced users

4.1 Curve analysis options

When the user clicks on a pixel in the map plotting area, a curve plot (intensity versus angle/image number) will be displayed in the second (right) plotting area. The scaling of the x axis of the curve can be selected via the "Delta-omega" input field in the top row of the visrock GUI: If a value of 1 is entered, the x axis can be interpreted as image numbers; if the angular increment used in the scan is entered, the x axis can be interpreted as angles. The "Start omega" input field can be used to select an angular offset (i.e. a shift of the entire x axis).

The exact meaning of the y -axis of the curve depends on the setting of the option selected in the curve analysis droplist just below the plot. The following functions are available:

orig data: Plot the intensity data contained at the respective pixel position in the original images, as a function of image number/angular position.

ROI-integral: Plot the sum of intensity values in all pixels visible in the map plotting area. When the map is in the completely unzoomed state, the sum of intensities in the entire image(s) is plotted.

ROI-min: Plot the minimum intensity value reached in any pixel across the current zoom area, as a function of image number/angular position.

ROI-max: Plot the maximum intensity value reached in any pixel across the current zoom area, as a function of image number/angular position.

spatial stddev: Plot the standard deviation of intensity values in a local surrounding of the pixel selected by clicking, as a function of image number. Useful e.g. for investigating whether the "noisyness" of a camera image changes with time or temperature. When this option is selected, an additional "stddev radius" input field will appear below the plot, allowing to select the width of the local surrounding (with a standard value of 1 pixel to the right and left).

sum_x(-log(data)) for ROI: Expert feature, undocumented so far.

sum_y(-log(data)) for ROI: Expert feature, undocumented so far.

4.2 Curve plotting options

The appearance of the Intensity-versus-angle curves displayed in the right drawing area can be changed via the buttons below the plot.




Figure 3: Curve plotting options

4.2.1 Plot range selection

The button labeled "log" allows to change back and forth between linear and logarithmic intensity scale.

The button labeled "Full range" determines the range of values shown on the y (intensity) axis: When the button is not pressed, the y range is individually adapted to the actual curve, thus maximizing the visible "contrast". When the button is pressed, the y range is determined from the minimum and maximum values in the entire 3D data block. This facilitates the comparison (on an absolute scale) of several different extracted curves, obtained by clicking in sequence on different pixels in the map.

4.2.2 Peak detection

The button "show peaks", when pressed, will activate automatic peak detection in the current curve. "Real" peaks (as opposed to local absolute maxima) will be marked by a sign. This serves mainly as a control mechanism for the "N_peaks" analysis option in map display.

Peak detection is based on a curve sharpening algorithm (see [Sil96]) followed by a search for local maxima in the sharpened curve. Maxima must be higher than adjacent minima by a factor that can be regulated via the "peak_detect.max_ratio" input field. The sharpened curve as an intermediate step towards peak detection can additionally be displayed by pressing the "show sharp" button.

4.2.3 Curve shape analysis

Curve shapes can be analyzed either by a phenomenological determination of half-widths (button `FWHM`) or by fitting with a linear (button `Lin.Fit`) or Gaussian (button `Gauss-Fit`) model. In all three cases, the resulting parameters will be displayed in a corner of the plot.

4.2.4 Noise reduction

The `median filter` button will apply a one-dimensional median filter to the current extracted curve. This will smoothen local spikes and reduced noise, but also flatten "real" peaks and (slightly) reduce resolution along the x direction. The width of the filter window can be regulated via the `smooth-/median-width` input field.

The `Dezinger` button will try to reduce noise from extreme intensity values in individual positions ("zingers") by applying a statistical criterion [BTG99]. This is mainly meant as a control facility for the `Dezinger` option in map plotting, which works according to the same algorithm.

4.2.5 Other plot options

Pressing the `oplot` button allows to plot a new curve without erasing the old one. This facilitates comparison of curves extracted from several different pixels.

When the `VBar` button is activated, a vertical red line in the curve plot indicates the angular position (image number in the series) at which the picture presently displayed in the map plot area was taken.

4.3 Image Plotting options

4.4 Axis units

The x and y axes in the map display can be labeled by physical units rather than pixel numbers. To this end, select `Edit plot options` from the Extra-functions menu. In the extra window that appears (see Fig. 4), select appropriate values for the variables `X_UNIT_NAME` and `X_SCALE_FACTOR` (the length of one pixel, in the unit defined by `X_UNIT_NAME`), and similarly for y .

The same window also allows to select the character size to be used in the two displays and in the legends, and to switch on (1) or off (0) the display of various titles and subtitles. Finally, the nature and format of the map legend display can be influenced via this window.

4.4.1 Color scale / Plot range selection

The color scale in the maps can be adjusted so as to maximize contrast and highlight specific features. Three mechanisms are available: One can either enter lower and upper limit values manually in the fields labeled `Minimum value:` and `Maximum value:`. As a consequence, any pixel values beyond these limits will be represented with the same




Figure 4: The "Edit Plot Options" window, allowing to define axis units, among others.

color as the limits themselves. In this way, the entire range of available colors (contrast) is concentrated on this interval of values.

As an alternative, the values can be defined interactively by using the two vertical sliders to the left of the map plotting area.

The third, most "automatic" way is to define min and max values in terms of percentiles: After entering the required numbers in the fields labeled "Minimum percentile" and "Maximum percentile" (see Fig. 5) the program will automatically translate these into limits in terms of absolute pixel values. Percentile values must lie between 0 and 100. As an example, min and max percentiles of 5 and 90 mean that the lowest 5% and the highest 10% of pixel values will be "flattened" (mapped to the corresponding limiting color). This is particularly useful when a small number of outlier pixels ("zingers") would otherwise spoil the automatic selection of color scales.




Figure 5: Map plotting options (selection).

4.4.2 Linear/logarithmic scaling. Aspect ratio

When the button "Map: log" is pressed, the color scale of the map shown in the left plotting area is chosen on a logarithmic scale. This allows to obtain more image contrast in regions of weak intensity. Choosing log-scaling will cause problems when parts or all of the map are in the negative-value range (as may happen for position-of-maximum maps, for example).

When the button is not active, the color scale is chosen linearly; this will enhance the visibility of the strongest peaks.

The button just below it, labeled "1:1", will activate the "maintain-aspect-ratio" function. This means that images with different numbers of pixels along the x and y axes will be plotted in a rectangular area, with dimensions on the screen (roughly) proportional to the respective number of pixels.

When this button is not active, images will be plotted in a quadratic regions on screen, often making better use of the available display space.

4.4.3 Zooming and unzooming

By default, all of the available image information is represented in the display. Details can be magnified by choosing an area of interest (AOI). An AOI is defined interactively by clicking with the mouse: Place the mouse pointer over one of the corners of a rectangular AOI, press the middle mouse button (or, in case your mouse has only two buttons, the left and right buttons simultaneously), move the mouse pointer to the diagonally opposed corner of the AOI to be defined while keeping the button(s) pressed, and release. During the move, the present size of the AOI just being defined will be visible as a bright rectangle overlaid on the image. After releasing, the new (magnified) plot will automatically be displayed on screen.

The center of the zoomed area can be shifted in the horizontal, vertical or any diagonal direction by using the "transfer buttons" (see block of 3×3 buttons in the lower right of Fig. 3). This will leave the magnification (zoom factor) unchanged. An exception is the central button (labeled "o"), which will not shift, but increase magnification by a factor 2 in both directions.

The inverse operation – unzooming, i.e. demagnification by a factor of 2 – is performed by clicking once with the middle mouse button.¹ The clicked pixel will become the center

¹More exactly: Clicking and releasing on the same pixel, without intervening motion.

of the new image, which will be twice as wide and high as before (provided that enough data are available).

A double-click with the middle mouse button will perform a complete unzoom, i.e. all of the available data will be displayed on screen.




Figure 6: Arrow buttons for shifting the region of interest

In order to shift the region of interest to an adjacent area without changing the scale (magnification factor), use the "arrow" buttons (Fig. 6). The outer 8 buttons allow to shift the region horizontally, vertically or diagonally. The central button (labeled "o") will zoom in by a factor of 2.

4.4.4 Mirroring / Image inversion

The map display can be mirrored in either the x- or the y-direction by activating the "Invert x" and/or "Invert y" buttons (see Fig. 5).²

4.4.5 Image rotation

The map display can be rotated by an arbitrary angle, by entering a value in the field labeled "Rotate by (deg):" – see Fig. 5.

Note that in the case of rotation by "odd" angles parts of the display area may be void for lack of data, will parts of the available data disappear beyond the border of the plotting area.

4.4.6 2D Fourier transforms

Instead of the "original" data, their 2D Fourier transform data can be displayed in the map plotting area. Fourier data are complex numbers, so the user can choose to plot the real or imaginary part, or their absolute values. This choice is made via a droplist on the right of the map plotting options (see Fig. 5). **TODO: Fourier-droplist abbilden.**

²Note that when any one of these two buttons is active, value and profile extraction from the map by mouse clicking will not work (at least in the present version of visrock).

4.5 Histogram display




Figure 7: Histogram display

The "Show histogram" button (see Fig. 1) allows to display a statistic of pixel values and their frequency of occurrence in the current image (i.e. a histogram of pixel values in the plotting area). The histogram is taken only over the currently zoomed region, and thus changes with the zoom state (center and magnification factor). A separate window (see Fig. 7) containing the histogram plot will be opened, which can be kept open for later use if required since it does not block action in the main window. A few buttons allow to manipulate the histogram plot:

"Save this plot": Output the plot to a postscript file. A default filename of "tmp.ps" will be suggested, but can be changed before final confirmation.

"Save these data": Store the plot data to an ASCII data file. A default filename of "tmp.dat" will be suggested, but can be changed before final confirmation. The file will contain two columns of data, and one row for each data point in the plot.

"Real data": Plot the original histogram data (default).

"Fourier data": Plot the 1D Fourier transform of the original data (rarely useful for histogram data).

xlog: Logarithmic x-axis (axis of pixel values).

ylog: Logarithmic y-axis (axis of frequency of occurrence).




Figure 8: Some map analysis functions.

FWHM: Automatically analyze the width of the highest peak in the plot, and write the result to a corner of the plot

”min” and ”max”: Specify a restricted range for the y-axis, in absolute frequency-of-occurrence numbers.

4.6 Line profile extraction from images

One-dimensional value profiles can be extracted interactively from the 2D map plots via mouse clicks: A double-click with the left mouse button will extract a horizontal line profile (z values as a function of x coordinate) from the map, crossing the clicked pixel. A double-click with the right mouse button will extract a vertical line profile (z values as a function of y coordinate).

4.7 Correlations

Correlation coefficients measure the degree of similarity between two images. They can have values between 0 (no similarity) and 1 (complete correspondence).

Correlation *functions* are a more general concept, telling whether two images can be made to be similar after a lateral (2D) shift.

4.7.1 Correlation coefficients and R-values

The button labeled **”Correlation coefficients”** will trigger the calculation of correlation coefficients between any pairs of images in the series of N images currently loaded into memory. The result is a (real and symmetric) $N \times N$ matrix of correlation coefficients, which will be displayed as a 2D map in a separate window (see Fig. 9).

This function is useful e.g. when inspecting a series of flat field images in tomography, which should all be representing the identical distribution. Considerable **”breaks”** in the 2D field of correlation coefficients may indicate an uncontrolled change in experimental conditions during exposure of the flat field series.

A similar concept, well-known to crystallographers, are the so-called **”R-values”** between pairs of images. They are calculated by pixel-wise summation of relative differences between pixel values in the two images. One of the images is defined as a reference, while the other is the **”sample”**. Therefore, the resulting $N \times N$ matrix is asymmetric. Apart from this, the result is mostly equivalent to the case of correlation coefficients, with the difference that R-values are additionally sensitive to changes in overall intensity values.

This function may therefore be useful to detect e.g. a sudden partial beam loss during exposure of a flat field image series.

4.7.2 Correlation functions

A correlation coefficient between a pair of images is a single (scalar) number, whereas a correlation function for a pair of digital images is again a 2D distribution, with the same number of pixels in x and y as the original images. The central element has the same value (between 0 and 1) as the "trivial" correlation coefficient, while the value of all other pixels indicate how much similarity could be achieved by shifting one of the two images by the corresponding number of pixels with respect to the other.

Clicking on the button labeled "Correlation functions" allows to calculate a correlation function between the currently selected image and a second one, to be selected from a list of filenames presented in a new window (see Fig. 11). After confirmation and calculation (which may take some time), the resulting 2D correlation function is saved as a postscript file, which will be displayed on screen using `ghostview` (on Linux/UNIX only). As the correlation function has complex values, there will be several plots presenting the real and imaginary components as well as the absolute values as 2D distributions, and several 1D projections of the same data (see Fig. 12).

4.8 Noise reduction

Fig. 13 shows the various options available for suppressing noise in 2D images inside `visrock`.

4.8.1 2D Median filtering

Median filtering is a well-known technique for suppressing "extreme" pixel values in 2D images. Every pixel value is replaced by the median of all pixel values inside a square box surrounding the respective pixel. The width of this box can be selected via the field labeled "Median: x/y ", and must have a value of 2 or higher. The default value of 0 indicates that no median filtering will be applied.

Median filtering is rather efficient at eliminating spurious features (single outliers or small groups of outlier pixels). It has the additional advantage of having rather limited "collateral effects", as compared to other averaging and smoothing techniques. It will, however, slightly degrade the lateral resolution in the image, especially at larger box widths.

2D median filtering is a reversible technique inside `visrock` in the sense that it is applied "on-the-fly" only to the intermediate data just being plotted to screen (or saved to file), i.e. to the image analysis results, and not to the original image data residing in memory. Therefore, reverting to a width value of 0 will reproduce exactly the original, noisy data.

This convenience comes at a price in terms of efficiency: The median filtering operation is carried out anew every time that the image is being plotted. This takes a little time (on




Figure 9: Map of correlation coefficients between pairs of images from a series of 10 primary beam images. Both the x and the y axis are in terms of image numbers, defined from 0 to $N - 1$. Note that the distribution is symmetric, and that values on the diagonal are equal to 1.00 – simply indicating that each image is identical to itself. The non-diagonal values show that images nr. 0–7 are rather similar to each other, whereas images 8 and 9 deviate more strongly. The 2×2 fields in the upper right corner show that these last 2 images are again rather similar to each other. In summary, beam conditions stayed constant until image nr. 7, after which an (unknown) sudden slight change occurred. The nature of this change may be further elucidated by other functions of the `visrock` program. As an example, the correlation *function* between e.g. images 6 and 8 may give hints whether a (lateral or vertical) shift of beam or camera gave rise to the change.




Figure 10: Map of "spatial R-values" between pairs of images from a series of 12 dark field images. Note again that values on the diagonal are trivial – with an R-value of 0.00 being equivalent to a correlation coefficient of 1.00, both indicating identity. Note also that unlike the previous figure, this distribution is not symmetric.




Figure 11: File selection widget, allowing to select a second image with which the one presently displayed is to be correlated.

the order of seconds). To accelerate the procedure, the median filtering can alternatively be applied to the original image data by pressing the "Apply to all data" button after entering the appropriate box width values. Be aware that this operation is irreversible (unless you are willing to read all the images again from hard disk).

4.8.2 3D median filtering


When using "Apply to all data" for median filtering of the original image data, an additional and independent box width can be entered into the field labeled "z=". When a value ≥ 2 is entered, an additional 1D median filtering operation will be applied to the pixel value array along the z direction (the "rocking curve") at each 2D pixel location. This operation relates the values in adjacent images to each other, unlike the 2D filtering operation that is applied to each image separately (and is carried out first).

4.8.3 "Dezinger"


The field labeled "Dezinger", when filled with a value > 0 , will activate the "Dezinger" function. This is a noise reduction technique that tries to eliminate spurious high pixel values ("zingers") based on a standard deviation criterion. The value entered roughly represents a multiple of standard deviations that can be tolerated as "statistical" deviations; more highly deviating values will be eliminated.

For this to work, a set of N (at least 2) independent images theoretically representing the identical "true" image, but distinct due to noise, are required. This function is particularly useful for treating dark field or primary beam images.

The algorithm is based on [BTG99], assumes spatial homogeneity of "noisyness", and works iteratively until convergence is achieved and a best compromise image is found. The calculation works independently for each pixel in the images; adjacent pixels are not related to each other.



Correlation function of /home/luebbert/Data/1999/06_99/caf-C.0006.edf.gz vs. /home/luebbert/Data/1999/06_99/caf-C.0008.edf.gz vs. /home/luebbert/Data/1999/06_99/caf-C.0006.edf.gz vs. /home/luebbert/Data/1999/06_99/caf-C.0008.edf.gz



Correlation function of /home/luebbert/Data/1999/06_99/caf-C.0006.edf.gz vs. /home/luebbert/Data/1999/06_99/caf-C.0008.edf.gz

Figure 12: Result of calculation of correlation function between a pair of images: Full correlation function (real part), full correlation function (absolute value), and projection onto the horizontal axis.

4.8.4 3D variance filter

Unlike Dezingering, the "3D variance filter" function for noise reduction will combine the values of several spatially adjacent pixels: For each pixel, the average and the standard deviation of pixel values in a local environment (cube of 3 pixels side length) is calculated and compared to the central pixel. If the pixel value deviates from the average by more than and (interactively entered) tolerance value times the standard deviation, then it is replaced by the mean.

This algorithm is realized in two alternatives: The "parallel" version treats all pixels identically, and always takes the original pixel values when calculating the local environment values. The "serial" version immediately takes the new values into consideration when treating the local environments for the next pixel, as it works its way through the entire data array.

Both variants are rather efficient at eliminating "real" outliers (and only those). However, they are pretty demanding in terms of computing time and (particular) working memory: Saving of intermediate results will require a maximum of 6.5 times the memory space occupied by the original image data. As an example, if you are treating an image




Figure 13: Noise reduction options.

series with 100 MB of data, and you don't have at least 650 MB of RAM, better don't use this function.

4.8.5 3D smoothing

A more basic way of 3D noise reduction is realized by the "Smooth:" fields: 3 numbers can be entered here, which represent the lengths along the 3 coordinate axes of a box around each pixel. Values inside this box are averaged, and the average value assigned as a new value to the central pixel. The result of this function is similar to using the "3D variance filter" with a tolerance level of 0, but more efficient in terms of execution time and memory requirements.

4.8.6 3D rebinning

When the "rebin" option rather than the "smooth" option is selected, the same smoothing mechanism will result in a new data block with reduced dimensions, i.e. smaller numbers of pixels along the three axes (by factors corresponding to the values entered into the three fields).

Note that these last two functions are irreversible, unless you are willing to read the original images again from hard disk.

4.9 Image processing options - sharpening, smoothing and peak detection

visrock can perform a number of image processing options on the data to be displayed in the map plotting area:

Original data: No image processing operation.

Deriv_x: Show (discrete) derivative of image data along the horizontal axis, instead of the original image. In other words, replace pixel value $u_{i,j}$ by $\partial_x u = u_{i+1,j} - u_{i,j}$.

Deriv_y: Similarly for the vertical direction

Modulus of Gradient: Absolute value of the 2D vector constructed from the two derivatives, i.e. $\sqrt{(\partial_x u)^2 + (\partial_y u)^2}$. Good for detecting areas of fast changes in image contents.

Phase of Gradient: Direction of the same vector, in $[0, 2\pi]$. Often not very useful.

Laplace 1: Convolution of original image with a 3×3 kernel given by

$$\begin{pmatrix} 0 & -1 & 0 \\ -1 & 4 & -1 \\ 0 & -1 & 0 \end{pmatrix}.$$

Good for detecting (local) peaks in the image.

Laplace 2: Convolution of original image with a 3×3 kernel given by

$$\begin{pmatrix} -1 & -1 & -1 \\ -1 & 8 & -1 \\ -1 & -1 & -1 \end{pmatrix}.$$

Also good for detecting (local) peaks in the image.

Laplace 3: Convolution of original image with a 3×3 kernel given by

$$\begin{pmatrix} 1 & -2 & 1 \\ -2 & 4 & -2 \\ 1 & -2 & 1 \end{pmatrix}.$$

A third way of detecting (local) peaks in the image. For details on the relevant differences between the three, see e.g. [Pra91].

Peaks: Smooth+Laplace1: First smooth image by a median filter of a given width, then perform peak detection (as in “Laplace 1”). Good for distinguishing noisy and spurious local peaks (which will be levelled out by the smoothing operation) from “real” ones.

mean-local: Convolution of original image with a 3×3 kernel given by

$$\frac{1}{8} \begin{pmatrix} 1 & 1 & 1 \\ 1 & 0 & 1 \\ 1 & 1 & 1 \end{pmatrix}.$$

Replaces each pixel by the average value of its 8 nearest neighbours. Good for image smoothing, and as a preparation for the **signal-mean-local** operation.

signal-mean-local: Similar, but additionally subtracting the value of the central pixel itself. Good for detecting prominent pixels that are very different from their local surroundings.

stddev-local: For each pixel, calculate the standard deviation of values of the 8 nearest neighbours, again by convolution. Good for detecting areas of strong local changes, and as a preparation for **SNR-local**.

SNR-local: For each pixel, calculate the local signal-to-noise ratio by dividing `signal-mean-local` by `stddev-local`. Helps to identify “real” (local) peaks in noisy images, by showing which pixels exceed their surroundings by more than the local noise level.

This list is of course very incomplete, compared to the many miraculous procedures contained in books like [Jäh02]. A few additions, like high-pass and low-pass filters, will be inserted in the near future. Further suggestions from users are always welcome.

5 Using visrock, III: Special functions for specific applications

visrock contains a few functions which are individually optimized for very specific applications, and may well be ignored by general users.

5.1 Functions for Rocking Curve Imaging

5.1.1 Dark field correction

In rocking curve imaging (as in tomography, and presumably most other fields of X-ray imaging) it is important to subtract an image of "dark counts" – caused by camera dark current, readout noise, intentional offset etc. – from each one of the actual experimental images. This task can be automated by using the field labeled "Dark field:" in the upper left of the visrock GUI.

The name of an image file containing a dark field image can either be entered directly into this field via the keyboard, or selected interactively from a list of files on hard disk via the button on the right. The dark field filename must be entered *before* reading the actual (non-dark) image files. The dark field data will be read into memory and displayed on screen for immediate check. Upon completion of the reading process of the main list of images, the dark field data will then be automatically subtracted from each and every image in the series before further processing. This may take some time, depending on the computer system being used.

Instead of a single dark field file, the user can alternatively specify a set of N dark field files simultaneously (e.g. by holding the SHIFT and/or CTRL-buttons pressed while selecting files). These will be understood as several (noisy) instances of the same basic image, and will be merged into a single "denoised" dark field using the "Dezingering" algorithm described in section 4.8.3.

5.1.2 Masking sample areas

In Rocking Curve Imaging, the diffracting sample areas often cover only part of the camera's field of view. The non-diffracting parts appear very noisy, particularly in the position-of-maximum and FWHM maps (since the analysis is carried out on spurious peaks). To mask out these areas (i.e. plot them in the background colour), the "Mask wafer area" can be used.

In order for this to work, the "wafer_present_yes_no" analysis option must have been selected first. This will carry out an analysis of each local rocking curve and decide whether it corresponds to diffracting material, based on two criteria: maximum Intensity and FWHM. These two quantities must lie between respective minimum and maximum values, which can be entered interactively. As a result, a binary map (1 = diffracting, 0 = non-diffracting) of the sample area will be displayed.

The limiting values can optionally be changed and analysis repeated (e.g. by pressing the "replot" button) until the binary image matches expectations. When this is the case, these data can be used for embellishing other plots by pressing the "Mask wafer area" button and then selecting the corresponding analysis option.

5.1.3 Correction of maxpos maps for dispersion effects

The position-of-maximum maps in Rocking Curve Imaging only partly reflect the sample curvature, since they also contain instrumental effects. In order to correct for this, the "Do corr. resol." button can be used. Note that this button is available only when the "maxpos" analysis option is active.

After activating the correction option, the program will ask for input data describing the experimental circumstances (monochromator, sample, reflections, geometry...). After confirmation, visrock will calculate a correction map, display it as a postscript file, and subtract it from the raw maxpos map. The correction can be undone by pressing the "Don't corr. resol." button.

The calculation is based on the following formulae:
(To be completed).

TODO: Explain angular correction

TODO: Bild einbinden

Figure 14: Input window for parameters describing the experimental geometry.

5.2 Functions for Tomography

visrock does *not* perform 3D reconstruction from projections. This can be done by programs like HST. However, visrock can be used for preparing experimental projection images for reconstruction by HST. This is particularly useful when working at beamlines which are not (yet) optimized for tomography, where unexpected misbehaviour of optics, rotation axes and/or detectors may prevent successful reconstruction unless it is diagnosed and documented.

5.2.1 Flat field calibration

A flat field correction can be performed automatically, similarly to the case of dark field correction described above. The name of one or several primary beam images must be entered into the field labeled "Prim. beam:" *before* reading the main series of images. After reading and (possibly) dark field correction, the experimental images will then be *divided* by the primary beam data. Note that this entails a change of data type, which in many cases means increased memory requirements (e.g.: changing from unsigned short to float means a doubling of RAM space).

5.2.2 Correction for primary beam decay

At most synchrotron sources, the primary beam intensity slowly decays between injections. To correct for this decay in the series of projection images, proceed as follows:

Select an area in the images which is illuminated by the beam, but never "sees" the sample. This can most easily be verified by switching to `minimum Intensity` analysis mode. Zoom to exactly this area, and switch to `ROI-Integral` curve analysis option. If the curve displayed in the curve drawing area is roughly constant, no correction is necessary. Otherwise, use "Normalize to ROI as background" function (in the "extra functions" droplist) to correct the entire set of images for the visible decay. After this, the curve on display should be exactly horizontal.

5.2.3 Sinogram display

Sinograms can be displayed by using the `x- ω` image display mode (for vertical rotation axis), or `y- ω` -mode (in the case of a horizontal rotation axis). Starting from an individual projection image, click anywhere into the image with the left mouse button. This will have the effect of fixing the `y` position (for `x- ω` -mode, or the `x` position for `y- ω` -mode) at which the sinogram is to be taken. Then, switch to the appropriate image analysis mode. Note that in these modes, the image-number slider (at the bottom of the `visrock` GUI) will change values and function, now allowing to switch to different `x` (or `y`) positions / sinogram numbers.

5.2.4 Finding the location of the rotation axis

To find the location (in pixel units) of the rotation axis (horizontal or vertical) in a stack of tomographic projections, proceed as follows:

In sinogram display, identify one strongly absorbing feature. Take the average between the `x` (or `y`) positions at $\omega = 0^\circ$ and $\omega = 180^\circ$.

5.2.5 Zero padding

In case the illuminated area in the experimental projections is not symmetrical around the (newly found) rotation axis, the image area can be extended by filling with zeroes (in logarithmic absorption mode, or equivalently ones in transmission mode). **TODO: Explain how!**

5.2.6 Detecting axis wobble

TODO: Discuss the extra functions in curve analysis (expert features...).

5.3 Functions for Fine-Slicing

The Fine-Slicing part of `visrock` is in many ways similar to the Beamish program [LSP00], although it is based on `xds` rather than on `mosflm`).

The main idea of the method is to collect rocking curves for many crystallographic reflections simultaneously. This is done by recording diffraction patterns with the standard rotation method and a CCD camera, but using a very small step width.

Each of the many Bragg spots needs to be located and indexed. This is done by an initial indexing step using the `xds` program [Kab02]. The remaining task carried out by `visrock` is to extract the rocking curve for each reflection from the individual images, apply the Lorentz correction and an additional geometry correction, output the reflection profiles, plus a big data array containing HKL, FWHMs, and many more.

After some initial enthusiasm, the Fine Slicing project was abandoned because it was found that the method did not provide more information than a set of classical rocking curves, but wasted more beam time (given an appropriate multi-circle diffractometer), hard disk space and processing resources. The documentation will therefore not be completed (although the corresponding functions in `visrock` do work).

5.4 Functions for Small-angle scattering and general Crystallography

For scattering patterns recorded either by the rotation method (in crystallographic data collection for structure determination) or by small-angle x-ray scattering (SAXS), it is often interesting to overlay the experimental images with *resolution rings*. These are circles around the primary beam spot, labeled with the value of their associated reciprocal-space resolution (in \AA^{-1}).

To activate the display of resolution rings in `visrock`, simply press the button labeled "show resolution rings". In order to obtain reasonable results, you should additionally define the appropriate scattering parameters by using the function "Edit resolution ring parameters" from the "Extra-functions" menu. This will open up an extra window (see Fig. 15) allowing to define the following values:

`ORIG_X / ORIG_Y`: Location of the primary beam spot in the image, in pixel units.

`WAVELENGTH`: X-ray wavelength used for recording the scattering pattern (in \AA).

`DISTANCE_MM`: Distance from sample to camera (in mm).

`PIXEL_SIZE_MM`: Size of one detector pixel (in mm).

5.5 Functions for Detector Calibration

Linear fitting - residual (dark field), slope (sensitivity), CHI^2 (linearity), `stddev` (reproducibility/noise)

Various types of X-ray cameras differ from each other not only in terms of pixel size, data format and readout speed, but also in terms of dark current, sensitivity, uniformity, linearity, and noise. The latter group of parameters can be determined experimentally from




Figure 15: Defining parameters for resolution ring display.

a series of pictures recorded with the camera in question, either under constant conditions or with variable X-ray intensity, exposure time, temperature etc.

This is a software manual, not a camera characterization handbook. We will therefore restrict the discussion to a few hints on how **visrock** can help in the analysis of such data, assuming that the appropriate measurements have already been performed.

The noise level in a series of theoretically identical images (e.g. a set of primary beam images recorded under identical conditions) can be evaluated by using the "standard deviation" image analysis option (see also p. 13). Probably even more useful is the option "noise \Leftrightarrow Poissonian?", which compares the experimental variance in the data to the one expected when assuming Poissonian noise. The resulting map will show where the variance is "real", i.e. caused by the incident beam statistics, and where it is "artificial", i.e. caused by deficient/noisy camera pixels.

Camera sensitivity, linearity and dark current can all be quantified from a series of

primary beam images recorded with increasing exposure time (or decreasing attenuation factor). The analysis of such data is then performed with the "Linear fit: ..." image analysis options. The following result maps are of interest:

Linear fit: Slope: This shows how the camera pixels react to increasing incident X-ray dose – i.e. the sensitivity, and its uniformity across the field of view.

Linear fit: Background: This option extrapolates the camera signal to zero incident intensity – i.e. it shows an alternative form of dark field, to be compared to the one recorded with closed shutter.

Linear fit: Chi²: This option shows whether or not the data set for each individual pixel could be fitted by a linear slope. If yes (low values of CHI^2), the respective pixel shows a linear reaction to increasing incident dose – which is most likely to be the desired behaviour.

6 Planned future developments of visrock

Here is an incomplete list of things that could, and maybe will, be implemented in **visrock** in the future:

- RC and extraction even for maps in inverted mode.
- More automatic normalization to monitor counts from *.edf file headers.
- Save memory space by cropping: Display first file of series, let user zoom to ROI, allocate space only for ROI
- 2D lattice tilt analysis (à la [MLK03]).
- Detector calibration: Check/correction for spatial distortion of tapered CCD cameras.
- Advanced 3D rendering of the entire data stack
- Multi-peak analysis (beyond the mere number of peaks per pixel).
- Segmentation of 3D data (identifying coherent crystallites, even if overlapping in space).

7 List of known bugs

7.1 Bugs

IDL crashes after you selected a large (≥ 200) number of individual files (via the "Read CCD files" button).

This appears to be an IDL bug, which occurs quasi predictably for file numbers above 250-300, at least on Linux systems. No solution has been found so far. A workaround is to select large series not manually, but via a text file containing the names of the image files (via the "Read file list" button).

Program stops while reading image series . The problem is currently under study. A workaround is not to touch the entire computer during reading (in particular, don't change between windows or virtual desktops).

Histogram display won't work, or percentile ranges are not translated into absolute display

In both cases, the histogram of map values cannot be calculated correctly by IDL. This problem is related to some IDL-internal functions. The reason for this failure is currently under study.

Solved with version of September 2004 - no longer a problem.

Program stops while producing an output file (postscript) Sometimes `visrock` stops working while writing a result image to postscript file, or (most often) while opening the postscript file for display on screen using `ghostview`. This is not really IDL crashing, but rather the `gv` child process that does not terminate. Unfortunately, the child process cannot be terminated externally, and IDL cannot be told to not wait for it either. The only solution: kill the entire IDL process, and start again.

To avoid this kind of crash, be a little patient: Take half a second between consecutive mouse clicks, and give the program time to "digest" your requests. This clearly helps to minimize problems of this kind.

7.2 Not bugs - features

In some cases please note that "It's not a bug – it's a feature". Known restrictions of this kind are:

- There is no way to stop a running movie – once you have started it, you must watch it till the end. If you want to save time, switch to linear (non-log) mode.
- Line extraction from maps does not work when the "invert x/y" buttons are active.
- Line extraction does not work when using physical length scales for the x/y axes, rather than pixel numbers.

- It is not possible to add new images to a series already in memory. Rather, one has to read the entire series again.
- There is no way to output a series of images as a movie (*.mpg or *.avi format, for example). Although this is theoretically very easy in IDL, it requires a special license.

An alternative: Movie creation is very simple on Linux using the tools from the ImageMagick suite (`convert`), in conjunction with `zing` for initial conversion from *.edf to other formats.

Generally, please be aware that this program was written by a single scientist, not a group of professional programmers. The focus is on some advanced data processing functions and innovative scientific applications, not on extreme user-friendliness, routine application and ultimate stability.

Therefore, it may happen that the program behave in unforeseen ways, or even crashes. Watch the IDL output window (on Linux, this is identical to the terminal window from which you started `visrock`). Some program errors cause error messages to be output here – but some others don't. If you observe an error message, please copy-and-paste it into an email to the author, also describing as much as possible of the circumstances. If you get an error, but no message – then be patient and try again. If (and when) you are convinced that it is a reproducible bug, please send a short email to the author so it can added to the list of known bugs.

8 Impressum

Author: Daniel Lübbert, Daniel.Luebbert@gmx.de

Send your comments and remarks, questions (if any), suggestions for improvement and donations directly to the author.

Send your complaints to `/dev/null`³.

³If you are not a Linux guru, please note: This is a joke.

9 Literature references

References

- [BSL00] H. D. BELLAMY, E. H. SNELL, J. LOVELACE, M. POKROSS, and G. E. O. BORGSTAHL. *The high-mosaicity illusion: revealing the true physical characteristics of macromolecular crystals*. Acta Cryst. D **56**(8), p. 986–995 (2000).
- [BTG99] S. L. BARNA, M. W. TATE, S. M. GRUNER, and E. F. EIKENBERRY. *Calibration procedures for charge-coupled device x-ray detectors*. Rev. Sci. Instr. **70**(7), p. 2927–2934 (1999).
- [Her80] G. HERMAN. *Image reconstruction from projections*. Academic Press New York (1980).
- [Jäh02] B. JÄHNE. *Digital Image Processing*. Springer Heidelberg 5th edition (2002).
- [Kab02] W. KABSCH. *XDS Program package* (2002). Available online at <http://www.mpimf-heidelberg.mpg.de/~kabsch/xds/>.
- [LBH00] D. LÜBBERT, T. BAUMBACH, J. HÄRTWIG, E. BOLLER, and E. PERNOT. *μm-resolved high resolution X-ray diffraction imaging for semiconductor quality control*. Nucl. Instr. Meth. B **160**(4), p. 521–527 (2000).
- [LM96] A. R. LANG and A. P. W. MAKEPEACE. *Reticulography: a simple and sensitive technique for mapping misorientations in single crystals*. J.Synchrotron Rad. **3**(6), p. 313–315 (1996).
- [LSP00] J. LOVELACE, E. H. SNELL, M. POKROSS, A. S. ARVAI, C. NIELSEN, N.-H. XUONG, H. D. BELLAMY, and G. E. O. BORGSTAHL. *BEAM-ish: a graphical user interface for the physical characterization of macromolecular crystals*. J.Appl.Cryst. **33**(4), p. 1187–1188 (2000).
- [MLK03] P. MIKULÍK, D. LÜBBERT, D. KORYTÁR, P. PERNOT, and T. BAUMBACH. *Synchrotron area diffractometry as a tool for spatial high-resolution three-dimensional lattice misorientation mapping*. J.Phys.D:Appl.Phys. **36**(10), p. A74–A78 (2003).
- [Pag03] E. PAGOT. *Diffraction Enhanced Imaging applied to Breast Tissue (?)*. PhD thesis Université Joseph Fourier, Grenoble (2003).
- [Pra91] W. PRATT. *Digital Image Processing*. Wiley New York 2nd edition (1991).
- [Sil96] Z. SILAGADZE. *A new algorithm for photopeak searches*. Nucl. Instr. Meth. A **376**, p. 451–454 (1996).