

Long Wavelength Phasing at SLS

Meitian Wang

Swiss Light Source


Paul Scherrer Institute

Outline

- Swiss Light Source and PX beamlines
- X06DA: a Beamline for Phasing Experiment
- Better Data for Phasing
- Phosphor-SAD Phasing

Swiss Light Source

- 2.4 GeV third generation synchrotron
- Top-up mode - 400 mA (± 1 mA)
- Low-emittance - $\epsilon_x = 5.5$ nmrad
- 18 beamlines




Protein Crystallography Beamlines

PX beamline	X06SA	X10SA	X06DA
Source	U19	U19	2.9T Superbend
Energy range (keV)	6.0 – 17.5	5.5 – 20.0	6.0 – 17.5
X-ray Optics	Dynamic focusing with adaptive optics		Achromatic fixed focusing
Focused beam size (μm)	HR: 85 × 10 MD2: 25 × 5	50 × 10	90 × 45
Divergence (mrad)	~ 0.4 × 0.06	0.6 × 0.1	2.0 × 0.5
Flux at 12.4 keV (photons/second)	2×10^{12}	2×10^{12}	5×10^{11}
Detector	PILATUS 6M	PILATUS 6M	PILATUS 2M
Sample changer	IRELEC CATS	CATS	CATS


X06DA: A Beamline for Phasing

- Bright super-bending magnet source
- True fixed-exit monochromator
- Achromatic focusing
- PILATUS 2M FAST (Nov. 2011)
- Multi-axis goniometer (Jan. 2012)
- *In situ* diffraction screening with SBS plate
- Integrated crystallization facility



SuperBend 2.9T: Bright and Stable




Stable and Rapidly Tunable Optics



Double Channel-Cut Monochromator



X06DA End-Station



PILATUS 2M FAST (60 Hz)

Excellent Data


SUBSET OF INTENSITY DATA WITH SIGNAL/NOISE ≥ -3.0 AS FUNCTION OF RESOLUTION

RESOLUTION LIMIT	NUMBER OF REFLECTIONS			COMPLETENESS OF DATA	R-FACTOR observed	R-FACTOR expected	COMPARED	I/SIGMA	R-meas	Rmrgd-F	Anomal Corr	SigAno	Nano
	OBSERVED	UNIQUE	POSSIBLE										
5.01	3702	380	382	99.5%	1.1%	1.2%	3701	170.73	1.2%	0.4%	68%	1.699	275
3.55	5798	646	646	100.0%	1.0%	1.2%	5798	164.48	1.1%	0.4%	35%	1.111	535
2.89	8225	834	834	100.0%	1.3%	1.4%	8225	149.63	1.3%	0.5%	32%	1.163	739
2.51	9896	966	966	100.0%	1.6%	1.6%	9896	124.92	1.7%	0.7%	29%	1.120	870
2.24	11029	1100	1100	100.0%	1.9%	1.9%	11029	102.49	2.0%	1.0%	17%	1.068	977
2.05	10859	1187	1187	100.0%	2.4%	2.2%	10859	82.11	2.5%	1.3%	11%	1.022	1072
1.90	13102	1326	1326	100.0%	3.3%	3.2%	13102	61.08	3.5%	2.1%	7%	0.941	1231
1.77	14043	1382	1382	100.0%	4.5%	4.5%	14043	45.70	4.7%	2.8%	7%	0.862	1285
1.67	15441	1497	1497	100.0%	6.4%	6.6%	15441	32.86	6.7%	3.9%	5%	0.840	1405
1.59	15253	1572	1572	100.0%	8.7%	8.7%	15253	24.43	9.2%	5.7%	3%	0.830	1412
1.51	15304	1636	1636	100.0%	12.6%	12.6%	15304	17.27	13.3%	8.9%	1%	0.796	1525
1.45	16779	1722	1722	100.0%	16.5%	16.9%	16779	13.17	17.4%	11.2%	2%	0.790	1625
1.39	17288	1799	1801	99.9%	22.9%	23.8%	17288	9.49	24.2%	15.5%	1%	0.775	1704
1.34	16702	1852	1853	99.9%	31.3%	32.8%	16701	6.70	33.2%	22.2%	-3%	0.738	1691
1.29	12555	1932	1933	99.9%	37.7%	39.8%	12551	4.58	41.0%	34.7%	-1%	0.740	1402
1.25	8302	1954	1960	99.7%	45.5%	48.7%	8222	2.88	51.9%	51.4%	-4%	0.699	949
1.22	5796	2010	2072	97.0%	53.2%	58.5%	5436	1.86	63.9%	73.8%	-1%	0.696	540
1.18	3505	1755	2108	83.3%	69.6%	75.5%	2863	1.14	89.6%	118.7%	-1%	0.612	134
1.15	1558	1124	2181	51.5%	68.2%	76.4%	801	0.87	93.1%	138.2%	-12%	0.549	14
1.12	338	303	2211	13.7%	71.3%	73.1%	70	0.71	100.9%	148.4%	0%	0.000	0
total	205475	26977	30369	88.8%	2.3%	2.4%	203362	35.60	2.4%	7.0%	3%	0.871	19385


- Insulin
- 10 Hz fine φ -slicing - 0.1 deg/0.1 sec
- 900 images, 90 degree data

Multi-Axis Goniometer (PRIGo)


- PRIGo: Parallel Robotics Inspired Goniometer
- Chi 0-90°, Phi 0-360°




Sphere of Confusion Measurement



PRIGo: Sphere of Confusion $< 5 \mu\text{m}$



Omega




Chi

PRIGo User Application




- Select best part of crystal
- Align long-cell axis to avoid overlap
- Align even-fold rotation axis to record Bijvoet pairs on same image
- Multi-pass data collection with different crystal orientation for true redundancy and better absorption correction

PRIGo: Select Best Part of Crystal



Move
Chi to 30 deg

Move to
crystal's edge








720 deg data collection with Chi = 30 deg




Session Edit View Bookmarks Settings Help

SUBSET OF INTENSITY DATA WITH SIGNAL/NOISE >= -3.0 AS FUNCTION OF RESOLUTION

RESOLUTION LIMIT	NUMBER OF REFLECTIONS OBSERVED	NUMBER OF REFLECTIONS UNIQUE	NUMBER OF REFLECTIONS POSSIBLE	COMPLETENESS OF DATA	R-FACTOR observed	R-FACTOR COMPARED expected	I/SIGMA	R-meas	Rmrgd-F	Anomal Corr	SigAno	Nano
5.41	18801	1287	1289	99.8%	2.3%	2.3%	18801	105.41	2.4%	0.8%	98%	495
3.84	30150	2278	2278	100.0%	2.4%	2.4%	30150	97.91	2.5%	0.9%	95%	991
3.14	37316	2962	2962	100.0%	3.0%	2.9%	37316	76.12	3.2%	1.3%	89%	1331
2.72	47827	3496	3496	100.0%	4.6%	4.4%	47827	53.66	4.8%	2.0%	84%	1599
2.44	55125	3924	3924	100.0%	6.1%							
2.22	59559	4367	4367	100.0%	8.4%							
2.06	57249	4768	4768	100.0%	11.0%							
1.93	65351	5078	5078	100.0%	18.2%							
1.82	69377	5294	5486	97.9%	27.7%							
total	448755	33454	33568	99.7%	5.0							




Jan. 27th 2012




Courtesy of E. Lorentzen, J. Basquin, MPG

PRIGo: Align Long-Cell Axis




- Eukaryotic ribosome 60S subunit
- Courtesy of S. Klinge, N. Ban, ETH



PRIGo: Bijvoet Pairs on Same Image



$P4_12_12$
47.41 47.41 70.44
Crystal oriented along a axis




In situ Diffraction Screening



- In routine user operation
- 2 min switch time
- Same centering and data collection interface


Easy Long-Wavelength Case

Sm^{3+} L_{III} -edge, 1.846 Å



Courtesy of N. Stefan, S. Streikov, Belgium

Challenge in S-SAD and P-SAD Phasing



- There is really no much anomalous signal!
- $\langle \Delta F^\pm \rangle / \langle F \rangle \approx 1\text{-}2\%$ for S-SAD at 1.6-2.0 Å
- $\langle \Delta F^\pm \rangle / \langle F \rangle \approx 2\text{-}3\%$ for P-SAD at 1.6-2.0 Å

Signal/Noise

$$\langle \Delta F^\pm \rangle / \langle \sigma(F) \rangle$$

$$\langle \Delta I^\pm \rangle / \langle \sigma(I) \rangle$$

$$I / \sigma(I)$$

Improve Signal to Noise

- Fine φ -slicing with PILATUS detector
- Dose fractionation with PILATUS detector
(Multiple data sets (MDS) strategy by B.C. Wang)
- Multi-pass data sets with different crystal orientation with PRIGo
- Sample mounting is as important as beamline instrumentation

Features of PILATUS Detector

- Point spread function of 1 pixel
- No dark current
- No readout noise
- Fast readout (2.3 ms) => shutter-free continuous rotation
- High frame rates (25 Hz and 60Hz)
- High dynamic range of 20 bits



Hybrid pixel: noise free detection

100 ms




- No readout-noise
- No dark-current.

1 hour




100K detector

Fine φ -slicing and background



Wide φ -slicing

- Large $\Delta\varphi$ ($\Delta\varphi > \xi$)
- Large overlap of reflections and background along φ
- Few images



Fine φ -slicing

- Small $\Delta\varphi$ ($\Delta\varphi \ll \xi$)
- Minimal overlap of reflections and background along φ
- Many images

Fine φ -slicing and Readout Noise

Detector with readout noise (CCD)



Detector without readout noise (PILATUS)



take full advantage of fine φ -slicing

minimize background

more images and smaller $\Delta\varphi$

The diagram for a PILATUS detector shows a horizontal blue arrow pointing to the right, labeled "more images and smaller $\Delta\varphi$ ". Above the arrow, the text "take full advantage of fine φ -slicing" is written in orange, and "minimize background" is written in green below it.

Overall statistics: $I/\sigma(I)$ and R_{merge}



- Test crystals: insulin, thaumatin, lysozyme, thermolysin
- 16 series of data sets with increasing $\Delta\phi$, each with 5 to 7 data sets
- 94 data sets in total

Dose Fractionation

Single data set approach:

$$\sigma_{total}^2 = G[I_s + I_{bg} + (m/n)I_{bg}] + m(K/A)^2 I_s^2$$

Counting statistics

Multiple Data set (MDS) approach:

$$\sigma_{total}^2 = G[I_s + I_{bg} + (m/n)I_{bg}] + \frac{m(K/A)^2 I_s^2}{N}$$









Systematic instrument error

This illustrates the theoretical advantage of applying the MDS approach over the traditional approach, which is summarized in the statement below:

For a fixed X-ray dose, collecting Multiple Data Sets with short exposures can produce better data than collecting a single set of data using the long exposures.

B.C Wang

Z.J. Lui, et al, Acta Cryst. (1011) A67, 544




A better data-collection strategy might have been *to increase the multiplicity while maintaining the total dose to avoid the inflation of the systematic error* term which limits the I/σ of individual observations. This would have resulted in higher I/σ values after averaging

Dose Fractionation with PILATUS

Detector with readout noise (CCD)








Detector without readout noise (PILATUS)




Better sample mounting

- - minimise plate bending



courtesy of J. Holton, ALS

Stress and Strain




courtesy of J. Holton, ALS

Two Examples of P-SAD

Hybrid RNA/ANA



- Hybrid decamer with one RNA strand and ANA (altritol nucleic acid)
- Unit cell $34.2 \times 34.2 \times 170.1 \text{ \AA}^3$ contains 1.5 hybrid decamer (15 base pairs \Rightarrow 30 phosphors) in asymmetric unit.
- *Sub-cell $34.2 \times 34.2 \times 34.0 (170/5) \text{ \AA}^3$ contains only 3 base pairs \Rightarrow 6 phosphors in asymmetric unit.*
- SHELXD with P-SAD sub-cell data set
- SHELXE with additional high-resolution data set (1.14 \AA), and -e0.914 to extrapolate to 0.914 \AA resolution

Hybrid RNA/ANA



Courtesy of M. Ovaere, L. van Meervelt, Belgium


Sarcin/Ricin Loop (SRL)



Sarcin/Ricin Loop (SRL), 26 P-site

SRL RNA Data

Estimated $\Delta F^\pm/F$	~ 2%
Wavelength (Å)	1.6
Space group	$P4_3$
Unit cell (Å)	29.59 x 29.59 x 76.09
Resolution (Å)	1.8 (1.9-1.8) \rightarrow 1.4
Oscillation range (°)	720
Completeness (%)	99.5 (97.3)
Redundancy	27.7 (23.4)
I / sigma	83.8 (52.2)
R_{meas}	0.041 (0.056)
$R_{\text{anom}} / R_{\text{p.i.m.}}$	0.014/0.008 = 1.75



SHELXD + SHARP + SOLOMON map


Sub-structure Determination

SRL RNA (26-site)	Crystal 1 P4 ₃ , 29.4 × 29.4 × 75.5 Å		Crystal 2
Rotation range	720 deg	360 deg	720 deg
Redundancy	28	14	28
R_{meas} (to 1.8 Å)	0.041 (0.056)	0.042 (0.061)	0.062 (0.163)
$R_{\text{anom}}/R_{\text{p.i.m.}} > 1.5?$	0.014/0.008 = 1.8	0.015/0.011 = 1.4	0.016/0.012 = 1.3
SHELXD: 2.1 Å cutoff	Partial solution with 15 correct sites	?	?
SHARP	Sub-struct. completion 26 sites	-	-
SHELXD: 1.8 Å cutoff	4 solutions in 1000 trials	5 solutions in 10'000 trials	?
SHELXE	Interpretable map	Interpretable map	-
SHARP	Excellent map	Good map	-


Phasing with Known P-Sites

SRL RNA (26-site)	Crystal 1			Crystal 2
Oscillation	720 deg	360 deg	180 deg	720 deg
Redundancy	28	14	7	28
R_{meas} (to 1.8 Å)	0.041 (0.056)	0.042 (0.061)	0.043 (0.059)	0.062 (0.163)
$R_{\text{anom}}/R_{\text{p.i.m.}}$	0.014/0.008= 1.8	0.015/0.011= 1.4	0.021/0.016= 1.3 <u>0.014/0.016= 0.9</u>	0.016/0.012=1.3
Sub-structure	26 correct sites			
PHASER+DM	Excellent map	Excellent map	Interpretable map	Excellent map
Map_CC	0.82	0.81	0.70	0.80

Phasing with Known P-Sites



X06DA for Phasing



Acknowledgment

Clemens Schulze-Briese (DECTRIS)

Vincent Olieric

Sandro Waltersperger

Marco Salathe

Claude Pradervand

Roman Schneider

Guanya Peng

Eric Ennifar

Philippe Dumas

SLS MX Group and Machine Group

B.C. Wang

our users !!!



Thank You for Your Attention!